


Instructions for use of the pressurized aerosol inhaler (PAI) with built-in spacer (*jet system*)

What is the *jet system*?


The aerosol cartridge is a device that contains a medication dissolved in the form of gas. By pressing the cartridge and breathing in through the mouthpiece, the medication enters the lungs with the air you breathe in. It includes a low volume circular spacer that reduces the disadvantages of the direct effect of the drug. The outer case is plastic and can be cleaned with soap and water.




How do you use the *jet system*?

Either standing or sitting down, relax as much as you can and follow these instructions:

1. Remove the lid.



2. Shake the inhaler vigorously.



3. Breathe out fully.

4. Lift your head and close your lips around the mouthpiece. Start breathing in through your mouth.




5. Press the inhaler and continue taking air in through your mouth, slowly and deeply (if you can't, breathe normally through the mouthpiece).

6. Hold your breath in your lungs for around 10 seconds, or for as long as you can.

7. Breathe out slowly.

8. Gargle and rinse your mouth after each use.

9. Replace the lid.



If you have any questions, please check with the nurse.